


UIC IST Teststand User Manual (v_4)

Yaping Wang
Zhenyu Ye

12/16/2012

1. IST Teststand




Note:

1) Local DAQ system has been installed at UIC for testing, and the cosmic trigger system is available for assembly.

2) The DAQ PC also provides service as a slow control server (slow control software - EPICS - should be installed, SL 5.3 installed).

3) The Iseg HV module is controlled through MPOD controller. The MPOD controller can be accessed via IP by EPICS software on DAQ server.


2. IST readout system



IST readout is identical to FGT design

- 1) For full readout testing of one IST stave, we need one ARM Readout Controller (ARC board) and two APV Readout Modules (ARM boards), and both of them are inserted in the VME crates.
- 2) FGT ARM Interface card (ABC board) is mounted on the backside of ARM board. The ABC board is used for inter-connection between Patch Panel Board (PPB) and ARM board. One ABC board covers two PPBs (connected by power cables and signal cables).
- 3) One IST stave has 3×12 APV chips, which were divided into three groups (A, B, C). They need one Transit Board (TB) and three PPBs for full readout. PPB and TB were connected with 50-channel bus.
- 4) The Data ReadOut Receiver Card (D-RORC, PCI kind) is mounted on the PCI slot of DAQ server, which builds a detector data link with the SIU card on the ARC via short optical fiber.
- 5) Slow control: (1) EPICS was installed on a Linux server (SL 5.3 suggested); (2) An EPICS package is available on the DAQ PC to control/monitor Iseg HV bias supply. Nominal bias voltage for sensors is set to 60 Volts.

3. IST VME crate




- 1) The ARC board should be located on the 7th slot of crate (If two ARC boards, the second one should be inserted into the 14th slot), and the SIU card is mounted on the ARC board.
- 2) The SIU card is connected to the D-RORC card via optical fiber.
- 3) Cosmic trigger is transmitted to the ARC trigger input for data acquisition.
- 4) The ARM board 0 locates on the 1st slot, then ARM1, ARM2, etc. Each ARC board covers 6 ARM boards.
- 5) The Iseg power supply locates on left slot close to MPODc module. The bias HV is set to 60 Volts.
- 6) The MPODc module should have a fixed IP address, which belongs to the same local network as the DAQ PC. We fixed a IP address for the MPODc, and also distributed a fixed IP for DAQ PC.

Note: Once the cards are available and cables are connected well, we can switch on the crate.

Do check carefully connections before you switch on the crate.

4. Cosmic trigger system



1) Cosmic trigger system consists of two sets of Scintillators + PMT, and the PMTs work at negative HV (The scintillator with bigger size is working at -1.820 kV, and the other smaller one is working at -1.48 kV).

Note: (1) Do not directly switch off the HV module once it works with high voltage outputs. Firstly, you should decrease the output voltage step by step. (2) Connections: the PMT HV input “Cathode” connects to the HV supply output; Connect the PMT output signal (BNC connector close to Cathode) to the Discriminator.

- (1) The HV for the scintillator 1 is supplied directly by channel 1 the High Voltage Regulated DC Power Supply (Model name HV-1547) at -1.82 kV.
- (2) The HV for the scintillator 2 is provided the combination of the Model HV-1547 and Model ES-7092 (Phototube High Voltage Zener Divider). Firstly, we input the -1.82 kV to the Model ES-7092 from output channel 2 of the model HV-1547, and then insert a pin into position 220 (black) at channel 2. Then the output voltage from the channel 2 of the Model ES-7092 is -1.60 kV which is supplied for the scintillator 2.
- (3) Check the HV connection carefully, and switch on power button if no question.
- (4) Adjust the voltage selection ($3 * 500V/step + 3 * 100V/step + 2 * 10V/step$).

2) The pulses output from PMTs are transmitted to LRS model 621BL (Discriminator).

3) The two readout signals from discriminators are sent to LRS model 365AL (DUAL FOUR-FOLD) for coincidence (AND logic calculation).

4) The coincidence output then sent to P\S model 794 (QUAD GATER) for trigger pulse width adjustment. **The trigger signal sent to ARC board should be 1 ms of pulse width.**

5. IP address setting for MPODc (Windows OS)

1) Get file or CD: "MPOD_CD_ROM_002-6"

2) Open the CD and go to sub-directory "NET SNMP", and install "net-snmpp-5.5.0-2.x64.exe". No special requirements for installation.

3) Once the net-snmpp installed successfully, go to sub-directory "MPOD_Software/MUSEcontrol", and install MUSEcontrolInstall-2.1.1563.0.exe

4) Start program "Muse control" and open menu "System"-->"Configure": Setting IP address, subnet mask. Using the default port number. **MPODc IP address is set to: 131.193.191.168**


5) Save and Exit the program.

6) Once the IP setting successfully for MPODc, we can check now: Open a web browser and input the MPODc's IP address, then Enter. We can monitor the status of Isege power supply's outputs.


Channel	Voltage	Current	Measured Sense Voltage	Measured Current	Measured Terminal Voltage	Status
U 0	0 V	10.000 mA	53.044 mV	0.0 uA	53.044 mV	OFF
U 1	0 V	10.000 mA	65.88 mV	0.0 uA	65.88 mV	OFF
U 2	0 V	10.000 mA	29.328 mV	0.0 uA	29.328 mV	OFF
U 3	0 V	10.000 mA	69.90 mV	0.0 uA	69.90 mV	OFF
U 4	0 V	10.000 mA	61.349 mV	0.0 uA	61.349 mV	OFF
U 5	0 V	10.000 mA	56.069 mV	0.0 uA	56.069 mV	OFF
U 6	0 V	10.000 mA	55.394 mV	0.0 uA	55.394 mV	OFF
U 7	0 V	10.000 mA	39.099 mV	0.0 uA	39.099 mV	OFF

6. Build the testing chain

1) Transition board (T-Board) and Patch Panel Board (PPB):



Backside of T-Board is mounted on the IST stave (A, B and C slots). Each slot on the T-Board front-side are connected to PPB board via a 50-channel bus. One T-Board covers to 3 PPBs.




The pin circled in yellow has to be connected to the central ground and as such also defines 0 for the bias voltage.


The two pins that are circled in blue are the bias connections. One pin is connected to the backplane connection of one sensor, the other pin is connected to the backplane of the other sensor.

2) ABC board

ABC board is an interface between PPB and ARM boards, which mounted on the backside of the ARM boards. Each ABC board take charges of communications with two PPB boards (two power cable + two sets of signal cables).



3) Test chain (Note: bias HV is set to 70 Volts in actual test)



6. EPICS software installation (SL linux OS)

1) Installation guide

0. log in starist.phy.uic.edu and run a terminal

```
su (to get root permission)
```

```
mkdir /opt/epics
```

```
cd /opt/epics
```

1. download base, extensions, and medm to current directory

```
http://www.aps.anl.gov/epics/download/index.php
```

2. setting variable:

```
export EPICS_HOST_ARCH=linux-x86_64
```

```
export EPICS=/opt/epics
```

3. cd /opt/epics

```
tar xvf baseR3.14.11.tar
```

```
ln -s base-3.14.11 base
```

```
tar xvf extensionsTop_20120904.tar
```

```
tar xvf medm3_1_7.tar
```

4. cd /opt/epics/extensions/src/

```
mv /opt/epics/medm3_1_7 ./medm
```

5. cd \$EPICS/base

```
make
```

6. yum install openmotif-devel.i686 openmotif-devel.x86_64

7. cd /opt/epics/extensions/src/medm

```
make
```

8. setting enviroment variables:

```
touch EPICS.sh and EPICS.csh
```

```
vim EPICS.sh
```

```
export EPICS_HOST_ARCH=linux-x86_64
```

```
export EPICS=/opt/epics
```

```
export PATH=$PATH:$EPICS/base/bin/linux-x86_64:$EPICS/extensions/bin/
```

```
linux-x86_64
```

```
vim EPICS.csh
```

```
setenv EPICS_HOST_ARCH linux-x86_64
```

```
setenv EPICS /opt/epics
```

```
setenv PATH $PATH:$EPICS/base/bin/linux-x86_64:$EPICS/extensions/bin/
```

```
linux-x86_64
```

9. sncseq-2.0.12 and asyn4-20 installation (download, extract, and make under the /opt/epics directory)

(1)

10. install net-snmp service

```
su (to get root permission)
```

yum install net-snmp (net-snmp-devel package is installed, then the commands “snmpd”, “net-snmp-config” and “snmptrapd” are available)

11. install net-snmp utils

download net-snmp-5.6.2.tar.gz to via link: <http://sourceforge.net/projects/net-snmp/files/net-snmp/5.6.2/>

```
tar zxvf net-snmp-5.6.2.tar.gz
```

```
cd net-snmp-5.6.2
```

```
./configure (default setting, no special requirement)
```

```
make
```

```
su
```

```
make install
```

12. test snmp and set community passwords

```
(1) snmpwalk -v 2c -m +WIENER-CRATE-MIB -c public 131.193.191.168  
SNMPv2-MIB::sysDescr.0 = STRING: WIENER MPOD (2688018, MPOD 2.1.1909.0,  
MPODslave 1.09, MPOD-BL 1.50 )
```

```
SNMPv2-MIB::sysObjectID.0 = OID: WIENER-CRATE-MIB::sysMainSwitch.0
```

```
DISMAN-EVENT-MIB::sysUpTimeInstance = Timeticks: (45460426) 5 days,  
6:16:44.26
```

```
SNMPv2-MIB::sysContact.0 = STRING:
```

```
SNMPv2-MIB::sysName.0 = STRING:
```

```
SNMPv2-MIB::sysLocation.0 = STRING:
```

```
SNMPv2-MIB::sysServices.0 = INTEGER: 79
```

```
(2) snmpwalk -v 2c -m +WIENER-CRATE-MIB -c guru 131.193.191.168  
snmpCommunityName
```

```
WIENER-CRATE-MIB::snmpCommunityName.public = STRING: "public"
```

```
WIENER-CRATE-MIB::snmpCommunityName.private = STRING: "private"
```

```
WIENER-CRATE-MIB::snmpCommunityName.admin = STRING: "admin"
```

```
WIENER-CRATE-MIB::snmpCommunityName.guru = STRING: "guru"
```

```
(3) snmpset -v 2c -m +WIENER-CRATE-MIB -c guru 131.193.191.168  
snmpCommunityName.guru s seCrET (setting of password)
```

```
(4) snmpwalk -v 2c -m +WIENER-CRATE-MIB -c seCrET 131.193.191.168  
snmpCommunityName
```

```
WIENER-CRATE-MIB::snmpCommunityName.public = STRING: "public"
```

```
WIENER-CRATE-MIB::snmpCommunityName.private = STRING: "private"
```

```
WIENER-CRATE-MIB::snmpCommunityName.admin = STRING: "admin"
```

```
WIENER-CRATE-MIB::snmpCommunityName.guru = STRING: "seCrET"
```

13. switch to your home directory

```
add environment variable: “setenv LD_LIBRARY_PATH /usr/local/lib:/usr/lib:  
$LD_LIBRARY_PATH” in .cshrc or “export LD_LIBRARY_PATH=  
$LD_LIBRARY_PATH:/usr/local/lib:/usr/lib” in .bashrc
```

```
add command line "source /opt/epics/EPICS.sh" in .bashrc or "source /opt/epics/  
EPICS.csh" in .cshrc
```

14. Install slow control package

```
(1) download the slow control package to the /home/ist/
```

- (2) mkdir ist_hv_devel
- (3) extract the package to the directory
- (4) cd ist_hv_devel
- (5) vim configure/RELEASE
 EPICS_BASE=/opt/epics/base
 TEMPLATE_TOP=\$(EPICS_BASE)/templates/makeBaseApp/top
 ASYN = /opt/epics/asyn4-20/
 SNCSEQ=/opt/epics/sncseq-2.0.12/
- (6) replace "fgttester/fgt_hv_devel" with "ist/ist_hv_devel" in all files
- (7) correct errors in snmpApp/src/Makefile
- (8) vim ist-test_display
 source setup_independent_env
 setenv EPICS_CA_AUTO_ADDR_LIST NO
 setenv EPICS_CA_ADDR_LIST "starist.phy.uic.edu"
 cd display
 medm -x FGTTestPowerSupply.adl &
- (9) vim setup_independent_env
 setenv EPICS_CA_SERVER_PORT 45582
 setenv EPICS_CA_AUTO_ADDR_LIST NO
 setenv EPICS_CA_ADDR_LIST "starist.phy.uic.edu"
- (10) mkdir /home/ist/.snmp
 mkdir /home/ist/.snmp/mibs/
 cp /home/ist/ist_hv_devel/WIENER-CRATE-MIB.txt /home/ist/.snmp/

mibs

- (11) cd iocBoot/iocsnmp
 vim ist-test.cmd
 (input Mpodc IP address and replace SYSTEM name "ist-test" with "isttest")
 vim ../../display/FGTTestPowerSupply.adl
 (replace SYSTEM name "ist-test" with "isttest")
 vim iocstart
 source ../../setup_independent_env
 ist-test.cmd

- (12) cd
 mkdir .snmp
 mkdir .snmp/mibs/
 cp /home/ist/ist_hv_devel/WIENER-CRATE-MIB.txt ~/.snmp/mibs/
 cp /home/ist/ist_hv_devel/ WIENER-CRATE-MIB.txt /usr/share/snmp/

mibs/

- (13) cd /home/ist/ist_hv_devel
 gmake clean
 gmake

Reference:

EPICS Installation guide referred links below:

<http://nibot-lab.livejournal.com/91787.html>

<https://pubweb.bnl.gov/~mdavidsaver/epics-doc/epics-starting.html>

7. Software preparation for DAQ and Analysis

IST DAQ PC: starist.phy.uic.edu

User account: ist

Password: XXXXXXXX

1) Software preparation for DAQ

1. log in IST DAQ PC
2. Start the afs service
3. open a terminal window
mkdir "cvs" to install software via afs service
cd
vim .bashrc (add following variables)
export PATH=\$PATH:
export CVSROOT=/afs/rhic.bnl.gov/star/packages/repository
export RTS_PROJECT=STAR
export PROJDIR=/RTS
export RTS_DEV_PREFIX=pro_
export USE_64BITS=1
4. change to "cvs" directory
cd cvs
klog XXXXX (AFS account at RACF)
cvs checkout StRoot/RTS
cvs checkout online/RTS

5. install driver for RORC64 (This step is already done by Tonko)

6. install DAQ system (compiled by Tonko)

```
cd cvs/online/RTS/src/FGT_TEST/  
cp /home/tonko/cvs/online/RTS/src/FGT_TEST/fgt_run_daq ./  
cp /home/tonko/cvs/online/RTS/src/FGT_TEST/fgt_run_daq.o ./  
cp /home/tonko/cvs/online/RTS/src/FGT_TEST/fgt_read_daq ./  
cp /home/tonko/cvs/online/RTS/src/FGT_TEST/fgt_read_daq.o ./  
cp /home/tonko/cvs/online/RTS/src/FGT_TEST/fgt_run_daq ./  
cp /home/tonko/cvs/online/RTS/src/FGT_TEST/ars_lib.o ./
```

2) Software preparation for Analysis

1. log in IST DAQ PC
2. In order to match the analysis library at STAR RCF, we need to make both the Scientific Linux version and gcc version identical to one at RCF. After that, we could share the identical library via AFS service (sl57_x8664_gcc432).
open a terminal window
mkdir gcc and download gcc-4.3.2 package (<http://gcc.gnu.org/mirrors.html>)
cd gcc/

```

tar zxvf gcc-4.3.2.tar.gz
mkdir gcc_build
cd gcc_build
su
../gcc-4.3.2/configure --prefix=/usr --enable-shared --enable-threads=posix --enable-
__cxa_atexit --enable-c99 --enable-long-long --enable-languages=c,c++ --disable-
libstdcxx-pch
make
make install
3. setting RCF environment
cd ~
vim .cshrc (add the following lines in to the script)
    setenv GROUP_DIR /afs/rhic.bnl.gov/rhstar/group
    if ( -r $GROUP_DIR/star_cshrc.csh ) then
        source $GROUP_DIR/star_cshrc.csh
    endif
4. open a terminal window and switch to your home directory
mkdir ARMdisplay
cd ARMdisplay
copy the latest script "ARMdisplay.C" to the current directory. You can find the latest
version at RCF: /star/institutions/mit/nieuwhzs/IST_TESTING/
IST_PROTOTYPE_STACK_02-03-04/COSMICS_Sep2012
Once the script is created first time or edited, you should compile it again:
rm ARMdisplay ARMdisplay.o
stardev
make
5. the executive file "ARMdisplay" is generated once compiling successfully.

```

8. Instruction to take data with IST

I. Starting readout crate control

- (1) Start EPICS software: the EPICS back-end software is running in the background on IST DAQ PC. **Once the IST DAQ PC gets rebooted, the EPICS can not be started automatically.**

login to the IST DAQ PC directly:

```
ssh ist@starist.phy.uic.edu
```

```
password: XXXXXXXX
```

```
cd /home/ist/ist_hv_devel/iocBoot/iocsnmp
```

```
./iocstart (let it run all the time)
```

- (2) Start the readout crate control GUI:


```
cd ~/ist_hv_devel
```

```
./ist-test_display
```


The GUI pops up now and parameters can be set easily.

- 1) Start the Main Switch (as shown in the following picture left)

click the switch, the circle turns green indicating the readout VME crate is powered on.



2) MpodC information can be found as shown in the picture below:



3) Enable the EHS 8205E HV output:

Click the button “V-EHS 8 06On”, and a new interface is popped up as below:



Firstly, we set the “Voltage” to 60.000

Secondly, switch on the v0 channel by click the switcher.

II. Running the DAQ

(1) login to the IST DAQ PC

go to the IST working directory

```
cd /cvs/online/RTS/src/FGT_TEST/
```

(2) when the readout crate is on and everything is hooked up correctly then you should be able to run the DAQ:

```
./fgt_run_daq -L -c ist_rdo_conf_YPWANG_for_test.txt -C 0x03 -n 100 -S /scratch/data/ist_software_test1.sfs
```

-L means that the ARC local clock is used, always needs to be there

-C 0x01 means that ARMO gets used (mask 0x01 for ARMO and 0x02 for ARM1, 0x03 for both ARMO and ARM1, etc.)
-c set a configuration file <configuration file>
-n set number of events
-S set output file name. The output should be put in /scratch/data

III. Analyzing DAQ files with ARMdisplay

- (1) login to IST DAQ PC
- (2) go to the ARMdisplay directory:
`cd ~/ARMdisplay/`
- (3) analyze daq file
`stardev`
`./ARMdisplay -f <daq_file_name>`
- (4) The ARMdisplay contains a lot of messages about missing APV's. It will also make four display windows and depending on how many ARM's were used up to 3 windows. The 4th display is for tracking results.
- (5) we can make changes for ARMdisplay.C, and "make" again.
Note: If more than one ARMs were used, we should make changes in the ARMdisplay.C, like readout slot, etc.
- (6) The cosmic data can be found here: /scratch/data/

